


Bacitracin Susceptibility Test



The **Bacitracin test** is used to distinguish **β -hemolytic streptococci**. It can distinguish between *Streptococcus pyogenes*, which forms a zone of inhibition around the bacitracin disc, and *Streptococcus agalactiae*, which grows up to the disc.

Results:

- Positive – *Streptococcus pyogenes*
- Negative – *Streptococcus agalactiae*

Links

Related articles

- CAMP test
- Catalase test
- Plasmacoagulase test
- Optochin test

Bibliography

- ŠMÍROVÁ, Václava. *Úvod do lékařské mikrobiologie*. - edition. -.
- BEDNÁŘ, Marek – SOUČEK, Andrej – FRÁŇKOVÁ, Věra, et al. *Lékařská mikrobiologie : Bakteriologie, virologie, parazitologie*. 1. edition. Praha : Marvil, 1996. 558 pp. ISBN 8023802976.

