

The Eames test



The **Eames test** is used to determine the plasticity period of amalgam. After mixing, the amalgam is divided into 4 parts and left to rest for 2.5 minutes. After 2.5 minutes, it is combined and mixed for 2 seconds, then divided again into 4 parts and left alone for 30 seconds. The procedure is repeated every 30 seconds until the moment when the individual pieces no longer come together during mixing.

Links

Related Articles

- Amalgam

References

- MAZÁNEK, Jiří – URBAN, František. *Stomatologické repetitorium*. 1. edition. Grada Publishing a.s, 2003. 456 pp. ISBN 80-7169-824-5.