

File:C-hmat.jpg



Size of this preview: 800 × 450 pixels.

Original file (2,304 × 1,296 pixels, file size: 279 KB, MIME type: image/jpeg)

How to hold an ambuvac mask

File history

Click on a date/time to view the file as it appeared at that time.

| | Date/Time | Thumbnail | Dimensions | User | Comment |
|---------|-------------------------------|-----------|------------------------|-------------------------------|---------|
| current | 00:22, 6 November 2022 | | 2,304 × 1,296 (279 KB) | JuneScndler (talk contribs) | |

You cannot overwrite this file.

File usage

The following page uses this file:

- Basic and advanced cardiopulmonary resuscitation

Metadata

This file contains additional information, probably added from the digital camera or scanner used to create or digitize it.

If the file has been modified from its original state, some details may not fully reflect the modified file.

| | |
|----------------------------------|-------------------------|
| Camera manufacturer | SONY |
| Camera model | HDR-XR160E |
| Exposure time | 1/25 sec (0.04) |
| F Number | f/2.4 |
| Date and time of data generation | 12:59, 14 November 2011 |
| Lens focal length | 8.2 mm |